

<b>Title:</b>	<b>Plate Sensor on Transfer Station</b>			<b>ID:</b>
				0169
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2004-12-02	2004-12-02	-	CMa	

## Q:

Got the carousel today – looks great. Just got one question – does the transfer station that you shipped to us for this system need a sensor attachment? We're just trying to work out if we need more holes in the hood base or not.

## A:

The transfer station plate detection sensor is optional. You may or may not use it. Usually it's a good idea to use it since it allows better plate tracing.