

<b>Title:</b>	<b>Manual Carrousel Rotation</b>			<b>ID:</b>
				0379
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2010-12-15	2010-12-15		C <i>Ma</i>	

**Q:**

Just to be sure... Is it OK if we manually rotate the carrousel to access the stackers while the unit is turned ON? What do you recommend?

**A:**

Sure, the incubators have a door switch that de-activates the carrousel on the open storages there is a switch to de-activate the carrousel.