

<b>Title:</b>	<b>System Busy Access Alarm</b>			<b>ID:</b>
				0394
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2014-08-30	2014-08-30	STX220	C <i>Ma</i>	

## Q:

I am qualifying a StoreX Incubator Model STX220. The manual indicates an alarm called "System Busy Access." What is this alarm and how can it be tested?

## A:

This is when you open the front door while the handler accesses plates - namely when the shovel is extracted.