

<b>Title:</b>	<b>Voltage Overveiw</b>			<b>ID:</b>
				125
<b>12.03.13</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
		-	CMa	

## ***Mains Voltage Overview***

Instrument			External Supply Voltage Range		Internal Voltage				Receptical(s)
Series	Capacity	Type	Voltage	Frequency	Climate		Handling		
LPX	44,110,220		100/110/208/230V	50/60Hz	-	-	85-264V	50/60Hz	6A
LPR	44,240		100/110/208/230V	50/60Hz	-	-	85-264V	50/60Hz	6A
LPT	220		100/110/208/230V	50/60Hz	-	-	85-264V	50/60Hz	6A
STX	44, 110, 220, 500	IC, HC, HR, DC	100/110/208/230V	50/60Hz	100/110V	60/50Hz	88-264V	47/63Hz	10A
	1000		100/110/208/230V	50/60Hz	100/110V	60/50Hz	88-264V	47/63Hz	10A/16A
	44, 110, 220, 500	DF	208/230V	50/60Hz	208/230V	60/50Hz	88-264V	47/63Hz	10A/16A
STR	44, 240	IC, HC, HR, DC	100/110/208/230V	50/60Hz	100/110V	60/50Hz	88-264V	47/63Hz	10A
	44, 240	DF	208/230V	50/60Hz	100/110V	60/50Hz	88-264V	47/63Hz	10A/16A
STT	602, 1000		208/230V	50/60Hz	203/230V	60/50Hz	88-264V	47/63Hz	10A/16A
	3k0		208/230V	50/60Hz	203/230V	60/50Hz	110/230V	60/50Hz	10A/16A
STC	all	ULT, DF	208/230V	50/60Hz	203/230V	60/50Hz	110/230V	60/50Hz	10A/16A
			208/400V	50/60Hz	208/400V	60/50Hz	110/230V	60/50Hz	10A/16A